

Accessories for Ready-to-use Minisarts and Re-usable Syringe Filter Holders



Dosing syringe

The dosing syringe is perfectly suitable in connection with a filter holder for a rapid and simple filtration. The new dosing syringe in combination with our Minisart filter holders are ideal for the wetting of our nutrient pad sets with sterile water.

The volume of the dosing syringe can be infinitely adjusted from 0.5 to 5.0 ml by turning the screw at the handle. The syringe is user- and maintenance-friendly. Moreover, it is very easy to handle and avoids so fatigue signs of the hand after longer use.

The dosing syringe can be disinfected by boiling out. It is not recommended to autoclave the syringes. If it is exceptionally necessary, the plastic handle has to be removed first.



3-way valve

Allows continuous filtration to be carried out, connected to a syringe and fitted on the outlet side with a filter holder.

Disposable syringes

They can be used with the 3-way valve and the filter holders with female luer lock inlet connection. One packet contains 12 individually packed needles and 12 disposable needles.

Needles

Fit on the luer slip outlets of the syringe filter holders. Allow silicone caps or rubber bungs to be pierced and the filtrate to be induced selectively into a tube or another vessel. The stainless steel needle is autoclavable.

Minicheck BP

Pressure gauge for the determination of the bubble point of Minisart syringe filter holders with luer lock connectors. For the determination of the bubble point of Sartolab P20 or Sartolab P20 plus units, the integrity test holder 18099 is required.

Order number for dosing syringe	16685--2	
Order number for 3-way valve	16639	Autoclavable (121°C).
Replacement parts	6986070	Sealing (4x)
	6986071	Pressure spring (2x)
	6986072	Fixing spring (2x)
	6986073	Perbunan valve (2x)
Order numbers for disposable syringes	16644E	5 ml volume, pack of 12
	16645E	10 ml volume, pack of 12
	16646E	20 ml volume, pack of 12
	16647E	50 ml volume, pack of 12
Order numbers for needles	01324	Stainless steel needle
	01325	Disposable needle
Order number for Minicheck BP	17799	